The eosAC-LT/LO automated flux chamber is purpose-built for applications that demand a large sampling area. The chamber design incorporates scientific best practices and features swappable components that allow for hassle-free repairs in the field. Seamless integration with leading gas analyzers and included flux analysis software enable simple and accurate measurements of trace gas flux.

- Field-Ready, Durable Design
- Large 52 cm Diameter Chamber Footprint
- Opaque (-LO) or Transparent (-LT) Options
- Field-Serviceable Components
- Seamless Integration with Leading Gas Analyzers
- Flux Analysis Software



Hardware		
Enclosure	Alum	inum/Acrylic/ABS
Gas Lines		Bev-A-Line® (PTFE optional)
Auxiliary Sensor Ports		2
Dimensions		
Chamber Volume (V)		0.072 m ³
Exposed Soil Area (SA)		0.21 m ²
Reach (chamber to analyzer)		30 m
Overall Size (max. L x W x H)		90 x 61 x 45 cm 90 x 61 x 92 cm
Weight		17 kg
Operating Environment		
Temperature		0 to 50 C
Humidity		0-100% RH non-condensing

NIO+O+ AII	anagifications	auhiaat	40	ahanaa	ith at	nation
Note. All	specifications	Subject	LO	change	WILLIOUL	HOUGE

Environmental Sensors				
Air/Chamber Temperature	0 to 50 C			
Barometric Pressure	30 to 110 kPa			
Auxiliary Sensors (Optional)				
Soil Temperature	-40 to 60 C			
Soil Moisture (VWC)	0.0 to 1.0 m ³ /m ³			
PAR	0 to 4000 $\mu mol\ m^{-2}\ s^{-1}$			
Power				
Operating Voltage	12-15 V DC			
Operating Power (in motion)	< 17 W (1.4 A)			
Operating Power (closed/open)	< 1 W (0.1 A)			
Other Features (Optional)				
Exterior	Transparent or Opaque			
Collar Height	10 cm 20 cm			
Base Extension	30 cm			